

Product datasheet for **TP726230**

NLRC3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NLR family CARD domain-containing protein 3/NLRC3
Species:	Human
Expression cDNA Clone or AA Sequence:	K641-Q852
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human NLR family CARD domain-containing protein 3/NLRC3is produced by E. coli. The target gene encodingK641-Q852 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	197358
UniProt ID:	Q7RTR2
Synonyms:	NLR family CARD domain-containing protein 3; NLRC3; CARD15-like protein; Caterpillar protein 16.2; NACHT, LRR and CARD domains-containing protein 3; Nucleotide-binding oligomerization domain protein 3
Protein Families:	Druggable Genome



[View online »](#)