

Product datasheet for **TP726238**

PPM1G Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein phosphatase 1G/PPM1G
Species:	Human
Expression cDNA Clone or AA Sequence:	G191-D546
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Protein phosphatase 1G/PPM1G is produced by E. coli. The target gene encoding G191-D546 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5496
UniProt ID:	O15355



[View online »](#)