

## Product datasheet for **TP726240**

### EHF Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ETS homologous factor/EHF
Species:	Human
Expression cDNA Clone or AA Sequence:	I2-T199
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human ETS homologous factor/EHFis produced by E. coli. The target gene encodingI2-T199 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	26298
UniProt ID:	<a href="#">Q9NZC4</a>
Synonyms:	ETS homologous factor; EHF; hEHF; ESE3; ESE3B; ESEJ; ETS domain-containing transcription factor; Epithelium-specific Ets transcription factor 3 (ESE-3)
Protein Families:	Transcription Factors



[View online »](#)