

Product datasheet for **TP726244**

Achaete scute homolog 3 (ASCL3) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Achaete-scute homolog 3/ASCL3
Species:	Human
Expression cDNA Clone or AA Sequence:	E41-S164
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Achaete-scute homolog 3/ASCL3is produced by E. coli. The target gene encodingE41-S164 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	56676
UniProt ID:	Q9NQ33
Synonyms:	Achaete-scute homolog 3; ASCL3; BHLHA42; ASH-3; hASH3; SGN1; Class A basic helix-loop-helix protein 42 (bHLHa42); bHLH transcriptional regulator Sgn-1



[View online »](#)