

Product datasheet for **TP726245**

Phospholipase D2 (PLD2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Phospholipase D2/PLD2
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-N141
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Phospholipase D2/PLD2is produced by E. coli. The target gene encodingT2-N141 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5338
UniProt ID:	O14939



[View online »](#)