

Product datasheet for **TP726247**

ERAS Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GTPase ERas/ERAS
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-P210
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human GTPase ERas/ERAS is produced by E. coli. The target gene encoding E2-P210 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	3266
UniProt ID:	Q7Z444
Synonyms:	GTPase ERas; ERAS; HRAS2; HRASP; Embryonic stem cell-expressed Ras
Protein Families:	Druggable Genome



[View online »](#)