

## Product datasheet for **TP726248**

### CCNO Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cyclin-O/CCNO
Species:	Human
Expression cDNA Clone or AA Sequence:	V2-L140
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Cyclin-O/CCNO is produced by E. coli. The target gene encoding V2-L140 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10309
UniProt ID:	<a href="#">P22674</a>
Synonyms:	Cyclin-O; CCNO



[View online »](#)