

## Product datasheet for **TP726249**

### **RPL34 Human Recombinant Protein**

#### **Product data:**

Product Type:	Recombinant Proteins
Description:	Recombinant Human Large ribosomal subunit protein eL34/RPL34
Species:	Human
Expression cDNA Clone or AA Sequence:	V2-K117
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Large ribosomal subunit protein eL34/RPL34is produced by E. coli. The target gene encodingV2-K117 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6164
UniProt ID:	<u><a href="#">P49207</a></u>



[View online »](#)