

Product datasheet for **TP726250**

DEPDC5 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GATOR complex protein DEPDC5
Species:	Human
Expression cDNA Clone or AA Sequence:	F380-A530
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human GATOR complex protein DEPDC5is produced by E. coli. The target gene encodingF380-A530 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	9681
UniProt ID:	O75140



[View online »](#)