

## Product datasheet for **TP726253**

### CTBP1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human C-terminal-binding protein 1/CTBP1
Species:	Human
Expression cDNA Clone or AA Sequence:	S320-L440
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human C-terminal-binding protein 1/CTBP1 is produced by E. coli. The target gene encoding S320-L440 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	1487
UniProt ID:	<a href="#">Q13363</a>



[View online »](#)