

Product datasheet for **TP726254**

Retinal protein 4 (UNC119) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein unc-119 homolog A/UNC119
Species:	Human
Expression cDNA Clone or AA Sequence:	K2-R160
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Protein unc-119 homolog A/UNC119is produced by E. coli. The target gene encodingK2-R160 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	9094
UniProt ID:	Q13432
Synonyms:	Protein unc-119 homolog A; Retinal protein 4 (hRG4); UNC119; RG4
Protein Families:	Druggable Genome, Stem cell - Pluripotency



[View online »](#)