

Product datasheet for **TP726258**

RNF11 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human RING finger protein 11/RNF11
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-N154
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human RING finger protein 11/RNF11is produced by E. coli. The target gene encodingG2-N154 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	26994
UniProt ID:	Q9Y3C5



[View online »](#)