

Product datasheet for **TP726259**

Synapsin III (SYN3) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Synapsin-3/SYN3 |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | G400-D580 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | Recombinant Human Synapsin-3/SYN3is produced by E. coli. The target gene encodingG400-D580 is expressed with a His tag at the C-terminus. |
| Storage: | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 8224 |
| UniProt ID: | Q14994 |
| Synonyms: | Synapsin-3; SYN3; Synapsin III |
| Protein Families: | Secreted Protein |



[View online »](#)