

Product datasheet for **TP726263**

KLF1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Krueppel-like factor 1/KLF1
Species:	Human
Expression cDNA Clone or AA Sequence:	G250-L362
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Krueppel-like factor 1/KLF1 is produced by E. coli. The target gene encoding G250-L362 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10661
UniProt ID:	Q13351
Synonyms:	Krueppel-like factor 1; Erythroid krueppel-like transcription factor (EKLF); KLF1; EKLF
Protein Families:	Transcription Factors



[View online »](#)