

Product datasheet for **TP726265**

GSH2 (GSX2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GS homeobox 2/GSX2
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-L304
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human GS homeobox 2/GSX2is produced by E. coli. The target gene encodingS2-L304 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	170825
UniProt ID:	Q9BZM3
Synonyms:	GS homeobox 2; Genetic-screened homeobox 2; Homeobox protein GSH-2; GSX2; GSH2



[View online »](#)