

Product datasheet for **TP726270**

RASGAP (RASA2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ras GTPase-activating protein 2/RASA2
Species:	Human
Expression cDNA Clone or AA Sequence:	Q750-S850
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Ras GTPase-activating protein 2/RASA2is produced by E. coli. The target gene encodingQ750-S850 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5922
UniProt ID:	Q15283
Synonyms:	Ras GTPase-activating protein 2; GTPase-activating protein 1m (GAP1m); RASA2; GAP1M; RASGAP
Protein Pathways:	MAPK signaling pathway



[View online »](#)