

Product datasheet for **TP726271**

ACTR3B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Actin-related protein 3B/ACTR3B
Species:	Human
Expression cDNA Clone or AA Sequence:	R156-T326
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Actin-related protein 3B/ACTR3Bis produced by E. coli. The target gene encodingR156-T326 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	57180
UniProt ID:	Q9P1U1
Synonyms:	Actin-related protein 3B; ACTR3B; ARP11; ARP4; ARP3-beta; Actin-like protein 3B; Actin-related protein ARP4



[View online »](#)