

## Product datasheet for **TP726272**

### ACSBG2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Long-chain-fatty-acid--CoA ligase ACSBG2/ACSBG2
Species:	Human
Expression cDNA Clone or AA Sequence:	L453-A618
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Long-chain-fatty-acid--CoA ligase ACSBG2/ACSBG2 is produced by E. coli. The target gene encodingL453-A618 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	81616
UniProt ID:	<a href="#">Q5FVE4</a>
Synonyms:	Long-chain-fatty-acid--CoA ligase ACSBG2; ACSBG2; BGR; UNQ2443/PRO5005; Acyl-CoA synthetase bubblegum family member 2; Arachidonate--CoA ligase ACSBG2; Bubblegum-related protein; PRTD-NY3



[View online »](#)