

Product datasheet for **TP726274**

RHOH Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Rho-related GTP-binding protein RhoH/RHOH
Species:	Human
Expression cDNA Clone or AA Sequence:	D48-F191
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Rho-related GTP-binding protein RhoH/RHOHis produced by E. coli. The target gene encodingD48-F191 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	399
UniProt ID:	Q15669
Synonyms:	Rho-related GTP-binding protein RhoH; RHOH; ARHH; TTF; GTP-binding protein TTF; Translocation three four protein
Protein Families:	Transcription Factors
Protein Pathways:	Leukocyte transendothelial migration



[View online »](#)