

Product datasheet for **TP726280**

APOL3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Apolipoprotein L3/APOL3
Species:	Human
Expression cDNA Clone or AA Sequence:	E75-I200
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Apolipoprotein L3/APOL3is produced by E. coli. The target gene encodingE75-I200 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	80833
UniProt ID:	O95236
Synonyms:	Apolipoprotein L3; APOL3; Apolipoprotein L-III (ApoL-III); TNF-inducible protein CG12-1 (CG12_1)



[View online »](#)