

Product datasheet for **TP726282**

LLGL2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human LLGL2 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	R620-E817
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human LLGL2 Protein is produced by E. coli. The target gene encoding R620-E817 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	3993
UniProt ID:	Q6P1M3
Synonyms:	LLGL scribble cell polarity complex component 2; HGL; Lethal(2) giant larvae protein homolog 2; LLGL2
Protein Families:	Druggable Genome
Protein Pathways:	Tight junction



[View online »](#)