

Product datasheet for **TP726283**

RASSF3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ras association domain-containing protein 3/RASSF3
Species:	Human
Expression cDNA Clone or AA Sequence:	K84-D238
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Ras association domain-containing protein 3/RASSF3 is produced by E. coli. The target gene encoding K84-D238 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	283349
UniProt ID:	Q86WH2
Synonyms:	Ras association domain-containing protein 3; RASSF3
Protein Families:	Druggable Genome



[View online »](#)