

## Product datasheet for **TP726284**

### ARRDC4 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Arrestin domain-containing protein 4/ARRDC4
Species:	Human
Expression cDNA Clone or AA Sequence:	S130-A230
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Arrestin domain-containing protein 4/ARRDC4 is produced by E. coli. The target gene encoding S130-A230 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	91947
UniProt ID:	<a href="#">Q8NCT1</a>
Synonyms:	Arrestin domain-containing protein 4; ARRDC4



[View online »](#)