

Product datasheet for **TP726285**

KLHL2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Kelch-like protein 2/KLHL2
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-P155
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Kelch-like protein 2/KLHL2is produced by E. coli. The target gene encodingE2-P155 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	11275
UniProt ID:	O95198
Synonyms:	Kelch-like protein 2; Actin-binding protein Mayven; KLHL2



[View online »](#)