

## Product datasheet for **TP726286**

### MTMR1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Myotubularin-related protein 1/MTMR1
Species:	Human
Expression cDNA Clone or AA Sequence:	D2-V110
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Myotubularin-related protein 1/MTMR1 is produced by E. coli. The target gene encoding D2-V110 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	8776
UniProt ID:	<u><a href="#">Q13613</a></u>



[View online »](#)