

Product datasheet for **TP726288**

OSBPL10 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human OSBPL10 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-A348
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human OSBPL10 Protein is produced by E. coli. The target gene encoding E2-A348 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	114884
UniProt ID:	Q9BXB5
Synonyms:	Oxysterol-binding protein-related protein 10; ORP-10; OSBP-related protein 10; OSBPL10; ORP10; OSBP9



[View online »](#)