

Product datasheet for **TP726289**

SCOC Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Short coiled-coil protein/SCOC
Species:	Human
Expression cDNA Clone or AA Sequence:	P45-K159
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Short coiled-coil protein/SCOC is produced by E. coli. The target gene encoding P45-K159 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	60592
UniProt ID:	Q9UIL1
Synonyms:	Short coiled-coil protein; SCOC; SCOCO; HRIHFB2072



[View online »](#)