

## Product datasheet for **TP726290**

### **PREI3 (MOB4) Human Recombinant Protein**

#### **Product data:**

|  |  |
|--|--|
| <b>Product Type:</b>                         | Recombinant Proteins   |
| <b>Description:</b>                          | Recombinant Human MOB-like protein phocein/MOB4  |
| <b>Species:</b>                              | Human  |
| <b>Expression cDNA Clone or AA Sequence:</b> | V2-A225  |
| <b>Tag:</b>                                  | C-His  |
| <b>Buffer:</b>                               | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4  |
| <b>Note:</b>                                 | Recombinant Human MOB-like protein phocein/MOB4is produced by E. coli. The target gene encodingV2-A225 is expressed with a His tag at the C-terminus.  |
| <b>Storage:</b>                              | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Stability:</b>                            | 12 months from date of despatch  |
| <b>Locus ID:</b>                             | 25843  |
| <b>UniProt ID:</b>                           | <a href="#">Q9Y3A3</a>   |
| <b>Synonyms:</b>                             | MOB-like protein phocein; MOB4; MOB3; MOBKL3; PHOCN; PREI3; CGI-95; 2C4D; Class II mMOB1; Mob1 homolog 3 (Mob3); Mps one binder kinase activator-like 3; Preimplantation protein 3                             |



[View online »](#)