

Product datasheet for **TP726291**

PGPEP1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Pyroglutamyl-peptidase 1/PGPEP1
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-H209
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Pyroglutamyl-peptidase 1/PGPEP1 is produced by E. coli. The target gene encoding E2-H209 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	54858
UniProt ID:	Q9NXJ5
Synonyms:	Pyroglutamyl-peptidase 1; PGPEP1; PGPI; 5-oxopropyl-peptidase; Pyroglutamyl aminopeptidase I (PAP-I); Pyroglutamyl-peptidase I (PGP-I); Pyrrolidone-carboxylate peptidase
Protein Families:	Druggable Genome, Protease



[View online »](#)