

Product datasheet for **TP726292**

MPP2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MAGUK p55 subfamily member 2/MPP2
Species:	Human
Expression cDNA Clone or AA Sequence:	E77-V254
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human MAGUK p55 subfamily member 2/MPP2is produced by E. coli. The target gene encodingE77-V254 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4355
UniProt ID:	Q14168
Synonyms:	MAGUK p55 subfamily member 2; MPP2; DLG2; Discs large homolog 2; Protein MPP2
Protein Families:	Druggable Genome



[View online »](#)