

Product datasheet for **TP726293**

SGSM3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Small G protein signaling modulator 3/SGSM3
Species:	Human
Expression cDNA Clone or AA Sequence:	F529-E724
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Small G protein signaling modulator 3/SGSM3 is produced by E. coli. The target gene encodingF529-E724 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	27352
UniProt ID:	Q96HU1
Synonyms:	Small G protein signaling modulator 3; SGSM3; MAP; RABGAPLP; RUTBC3; Merlin-associated protein; RUN and TBC1 domain-containing protein 3; Rab-GTPase-activating protein-like protein (RabGAPLP)



[View online »](#)