

## Product datasheet for **TP726299**

### FBXL20 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human F-box/LRR-repeat protein 20/FBXL20
Species:	Human
Expression cDNA Clone or AA Sequence:	R80-V280
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human F-box/LRR-repeat protein 20/FBXL20 is produced by E. coli. The target gene encodingR80-V280 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	84961
UniProt ID:	<a href="#">Q96IG2</a>
Synonyms:	F-box/LRR-repeat protein 20; F-box and leucine-rich repeat protein 20; F-box/LRR-repeat protein 2-like; FBXL20; FBL2
Protein Families:	Druggable Genome



[View online »](#)