

Product datasheet for **TP726300**

ARF4L (ARL4D) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ADP-ribosylation factor-like protein 4D/ARL4D
Species:	Human
Expression cDNA Clone or AA Sequence:	R92-R201
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human ADP-ribosylation factor-like protein 4D/ARL4Dis produced by E. coli. The target gene encodingR92-R201 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	379
UniProt ID:	P49703



[View online »](#)