

## Product datasheet for **TP726303**

### EXOC3 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Exocyst complex component 3/EXOC3
Species:	Human
Expression cDNA Clone or AA Sequence:	N459-D691
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Exocyst complex component 3/EXOC3is produced by E. coli. The target gene encodingN459-D691 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	11336
UniProt ID:	<u><a href="#">O60645</a></u>



[View online »](#)