

Product datasheet for **TP726304**

WDR73 Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human WD repeat-containing protein 73/WDR73 |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | L163-L366 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | Recombinant Human WD repeat-containing protein 73/WDR73is produced by E. coli. The target gene encodingL163-L366 is expressed with a His tag at the C-terminus. |
| Storage: | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 84942 |
| UniProt ID: | Q6P4I2 |
| Synonyms: | WD repeat-containing protein 73; WDR73; HSPC264 |



[View online »](#)