

Product datasheet for **TP726309**

Aconitase 1 (ACO1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cytoplasmic aconitate hydratase/ACO1
Species:	Human
Expression cDNA Clone or AA Sequence:	G110-D250
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Cytoplasmic aconitate hydratase/ACO1is produced by E. coli. The target gene encodingG110-D250 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	48
UniProt ID:	P21399



[View online »](#)