

Product datasheet for **TP726310**

BBOX1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Gamma-butyrobetaine dioxygenase/BBOX1
Species:	Human
Expression cDNA Clone or AA Sequence:	S1475-L1768
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Gamma-butyrobetaine dioxygenase/BBOX1 is produced by E. coli. The target gene encoding S1475-L1768 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	8424
UniProt ID:	O75936



[View online »](#)