

## Product datasheet for **TP726314**

### **DDB1 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human DNA damage-binding protein 1/DDB1
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	D271-L498
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	Recombinant Human DNA damage-binding protein 1/DDB1is produced by E. coli. The target gene encodingD271-L498 is expressed with a His tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	1642
<b>UniProt ID:</b>	<u><a href="#">Q16531</a></u>



[View online »](#)