

## Product datasheet for **TP726316**

### VPS35 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human VPS35 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-E250
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human VPS35 Protein is produced by E. coli. The target gene encoding A2-E250 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55737
UniProt ID:	<u><a href="#">Q96QK1</a></u>

[View online »](#)