

Product datasheet for **TP726317**

RRM1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human RRM1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	N310-C500
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human RRM1 Protein is produced by E. coli. The target gene encoding N310-C500 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6240
UniProt ID:	P23921



[View online »](#)