

## Product datasheet for **TP726321**

### **NPL4 (NPLOC4) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human NPLOC4 Protein
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	E2-G165
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	Recombinant Human NPLOC4 Protein is produced by E. coli. The target gene encoding E2-G165 is expressed with a His tag at the C-terminus.
<b>Storage:</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	55666
<b>UniProt ID:</b>	<a href="#">Q8TAT6</a>
<b>Synonyms:</b>	Nuclear protein localization protein 4 homolog; NPLOC4; Protein NPL4; NPLOC4; KIAA1499; NPL4



[View online »](#)