

Product datasheet for **TP726322**

ME1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NADP-dependent malic enzyme/ME1
Species:	Human
Expression cDNA Clone or AA Sequence:	R130-K250
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human NADP-dependent malic enzyme/ME1 is produced by E. coli. The target gene encodingR130-K250 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	4199
UniProt ID:	P48163



[View online »](#)