

Product datasheet for **TP726324**

ATP6V1B2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ATP6V1B2 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	D346-I538
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human ATP6V1B2 Protein is produced by E. coli. The target gene encoding D346-I538 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	526
UniProt ID:	P21281



[View online »](#)