

Product datasheet for **TP726330**

DPP3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Dipeptidyl peptidase 3/DPP3
Species:	Human
Expression cDNA Clone or AA Sequence:	D29-W223
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Dipeptidyl peptidase 3/DPP3 is produced by E. coli. The target gene encoding D29-W223 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10072
UniProt ID:	Q9NY33
Synonyms:	Dipeptidyl peptidase 3; DPP3; Dipeptidyl aminopeptidase III; Dipeptidyl arylamidase III; Dipeptidyl peptidase III (DPP III); Enkephalinase B
Protein Families:	Druggable Genome, Protease



[View online »](#)