

## Product datasheet for TP726331

### RAVER1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ribonucleoprotein PTB-binding 1/RAVER1
Species:	Human
Expression cDNA Clone or AA Sequence:	C2-P100
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Ribonucleoprotein PTB-binding 1/RAVER1 is produced by E. coli. The target gene encoding C2-P100 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	125950
UniProt ID:	<a href="#">Q8IY67</a>
Synonyms:	Ribonucleoprotein PTB-binding 1; RAVER1; KIAA1978; Protein raver-1



[View online »](#)