

Product datasheet for **TP726337**

SNRPD1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Small nuclear ribonucleoprotein Sm D1/SNRPD1
Species:	Human
Expression cDNA Clone or AA Sequence:	E92-Q239
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Small nuclear ribonucleoprotein Sm D1/SNRPD1is produced by E. coli. The target gene encodingE92-Q239 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	6632
UniProt ID:	P62314



[View online »](#)