

## Product datasheet for **TP726338**

### Plectin (PLEC) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Plectin/PLEC
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-K228
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Plectin/PLEC is produced by E. coli. The target gene encoding A2-K228 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5339
UniProt ID:	<a href="#">Q15149</a>
Synonyms:	Plectin; PLEC; PLEC1; PCN; PLTN; Hemidesmosomal protein 1 (HD1); Plectin-1
Protein Families:	Druggable Genome



[View online »](#)