

Product datasheet for **TP726340**

ANK1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ankyrin-1/ANK1
Species:	Human
Expression cDNA Clone or AA Sequence:	D589-F895
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Ankyrin-1/ANK1 is produced by E. coli. The target gene encodingD589-F895 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	286
UniProt ID:	P16157
Synonyms:	Ankyrin-1; ANK-1; ANK1; ANK; Ankyrin-R; Erythrocyte ankyrin
Protein Families:	Transmembrane



[View online »](#)