

Product datasheet for **TP726341**

TJP1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Tight junction protein ZO-1/TJP1
Species:	Human
Expression cDNA Clone or AA Sequence:	N439-N575
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Tight junction protein ZO-1/TJP1 is produced by E. coli. The target gene encoding N439-N575 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	7082
UniProt ID:	Q07157



[View online »](#)