

Product datasheet for **TP726342**

CLTCL1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Clathrin heavy chain 2/CLTCL1
Species:	Human
Expression cDNA Clone or AA Sequence:	V509-L858
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Clathrin heavy chain 2/CLTCL1 is produced by E. coli. The target gene encoding V509-L858 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	8218
UniProt ID:	P53675



[View online »](#)